

Product Datasheet

Human CD274 protein (orb420013)

Description	Recombinant of human CD274 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	C-His, variant
Protein Sequence	Phe19-Thr239
Source	Human Cells
Uniprot ID	Q9NZQ7
Expiration Date	6 months from date of receipt.