



Product Datasheet

ERBB3 antibody (orb420032)

Description	ERBB3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ERBB3. Please contact us for the exact immunogen sequence.
Target	ERBB3
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	148 kDa
Uniprot ID	P21860, Q61526, Q62799
NCBI	001005915
Entrez	2065
Dilution Range	IHC-P: 1:30-100, WB: 1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat lung (lane ...)



IHC-P image of mouse testis tissue using...