

Product Datasheet

SLC40A1 antibody (orb420040)

Description

SLC40A1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human SLC40A1. Please contact us for the exact immunogen sequence.

Target

SLC40A1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Immunohistochemical staining of human br...

Note

For research use only

Isotype


IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography



IHC-P staining of rat colon tissue using...

MW

63kDa

Uniprot ID
[Q923U9](#), [Q9JHI9](#), [Q9NP59](#)
NCBI
[014585](#), [055400](#)
Entrez
[30061](#)
Dilution Range

IHC-P: 1:100-1:600

Expiration Date

12 months from date of receipt.