



Product Datasheet

DDX58 antibody (orb420060)


Description	Goat polyclonal antibody to DDX58
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	DDX58
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Spleen, Placenta, and Brain (Cortex). Recommended concentration: 3.75µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and HeLa cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of U937 cells. Recommended concentration:
Protein Sequence	KDFHFEKIPFDPAE
Clonality	Polyclonal
MW	107
Entrez	23586
Dilution Range	ELISA: 1:32000, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Immunohistochemical staining of Human Pl...



Immunohistochemical staining of Human Co...