

Product Datasheet

MMP7 antibody (orb420072)

Description

MMP7 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunohistochemical staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human MMP7. Please contact us for the exact immunogen sequence.

Target

MMP7

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml & - 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW



30 kDa

Uniprot ID
P09237, Q10738
NCBI
002414, 002423
Entrez
4316
Dilution Range


IHC-P: 1:100-600, WB: 1:200-1000

Expiration Date

12 months from date of receipt.



IHC-P image of Guinea pig kidney tissue ...



IHC-P staining of Guinea pig kidney tiss...