

Product Datasheet

MYD88 antibody (orb420074)

Description

MYD88 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunohistochemical staining of Guinea p...

Immunogen

KLH conjugated synthetic peptide derived from human MYD88. Please contact us for the exact immunogen sequence.

Target

MYD88

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

IHC-P staining of mouse stomach tissue u...

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

Immunohistochemical staining of rat sple...

MW

33 kDa

Uniprot ID

[P22366](#), [Q6Y1S1](#), [Q99836](#)

NCBI

[001166037](#), [001172566](#)

Entrez

[4615](#)

Dilution Range

IHC-P: 1:100-300, WB: 1:100-500

Expiration Date

12 months from date of receipt.