



Product Datasheet

PEAK1 antibody (orb420075)


Description	PEAK1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PEAK1. Please contact us for the exact immunogen sequence.
Target	PEAK1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	193 kDa
Uniprot ID	D4A563 , Q9H792 , Q69Z38
NCBI	024776 , 079052
Entrez	79834
Dilution Range	IHC-P: 1:100-300, WB: 1:100-800
Expiration Date	12 months from date of receipt.



WB analysis of rat liver (lane 1), rat b...



IHC-P image of rat kidney tissue using P...



IHC-P staining of rat bone marrow tissue...