

Product Datasheet

ENPP1 antibody (orb420076)

Description

ENPP1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human ENPP1. Please contact us for the exact immunogen sequence.

Target

ENPP1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW


105 kDa

Uniprot ID
[Q924C3](#), [P22413](#), [P06802](#)
NCBI
[006208](#), [006199](#)
Entrez
[5167](#)
Dilution Range


IHC-P: 1:100-600, WB: 1:100-800

Expiration Date


12 months from date of receipt.



IHC-P image of rat spleen tissue using ENPP1 antibody.



IHC-P staining of rat kidney tissue using ENPP1 antibody.



Immunohistochemical staining of mouse liver tissue using ENPP1 antibody.