

Product Datasheet

Human Prolactin R / PRLR Protein (orb420091)

Description Human Prolactin R, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 25 - Asp 23...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Prolactin R

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.0 kDa. The protein migrates as 64 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity 95%


Source Human Prolactin R, Fc Tag (orb420091) is expressed from human 293 cells (HEK293). It contains AA Gln 25 - Asp 234 (Accession # P16471-1).

MW 51.0 kDa

Uniprot ID **P16471**

NCBI **NP_000940.1**

Expiration Date 6 months from date of receipt.



Immobilized
Human
Prolactin,
Mouse
IgG2a...