

Product Datasheet

Human TIGIT Protein (orb420096)

Description	T-cell immunoreceptor with Ig and ITIM domains (TIGIT) is also known as V-set and immunoglobulin dom...
Reactivity	Human
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	TIGIT
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-mFc(mIgG2a)
Note	For research use only
Application notes	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 40.1 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAI01290.1
Purity	95%
Source	Human TIGIT, Mouse IgG2a Fc Tag (orb420096) is expressed from human 293 cells (HEK293). It contains AA Met 22 - Pro 141 (Accession # Q495A1-1).
MW	40.1 kDa
NCBI	AAI01290.1
Expiration Date	6 months from date of receipt.