

Product Datasheet

1X PBS Buffer with 1% Casein (orb420117)

Description	1X PBS Buffer with 1% Casein
Conjugation	Unconjugated
Tested Applications	ELISA, FLISA, IP, WB
Preservatives	Preservative: None. Stabilizer: None. See application note.
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store container at 2° to 8° C prior to opening. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts.
Note	For research use only
Application notes	This product is a ready-to-use 1X solution and no further dilution is required.
Purity	1X PBS, pH 7.2 buffer with 1% casein was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Dilution Range	ELISA: User Optimized, FLISA: User Optimized, IP: User Optimized, WB: User Optimized
Expiration Date	12 months from date of receipt.