

Product Datasheet

1X TBS Buffer with 1% Casein (orb420119)

Description	1X TBS Buffer with 1% Casein
Conjugation	Unconjugated
Tested Applications	ELISA, FLISA, IP, WB
Preservatives	Proclin is added as an antimicrobial agent.
Form/Appearance	Liquid
Concentration	1 X
Storage	Store container at 2° to 8° C prior to opening. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts.
Note	For research use only
Application notes	This product is a ready-to-use 1X solution and no further dilution is required.
Purity	The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Dilution Range	ELISA: User Optimized, FLISA: User Optimized, IP: User Optimized, WB: User Optimized
Expiration Date	12 months from date of receipt.