

Product Datasheet

DNA Salmon Testes (orb420123)

Description	DNA Salmon Testes
Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None. None
Form/Appearance	Lyophilized
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Salmon Testes DNA is suitable for use as a substrate for deoxyribonuclease assays and other molecular biological techniques.
Purity	Prepared by a modification to the method of Emanuel, C.F., and Chaikoff, I.L., : JBC, 203, 164 (1953). ≥75% native nucleic acid. Supplied dried.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.