

Product Datasheet

DNA Lambda Phage (orb420125)

Description	DNA Lambda Phage
Conjugation	Unconjugated
Preservatives	None. See application note.
Form/Appearance	Liquid (sterile filtered)
Concentration	438 µg/ml
Storage	Store vial at -20° C or colder. Expiration date is six (6) months from date of opening vial.
Note	For research use only
Application notes	A solution in 10 mM Tris-Cl and 1 mM EDTA, pH 8.0. Note: A higher sample load may be required to clearly see the 564 and 125 base pair bands.)
Purity	Lambda DNA fragments prepared by the digestion of purified lambda DNA with the restriction endonuclease Hind III. On agarose gel electrophoresis the mixture separates into eight individual bands having the following number of base pairs: 23130, 9416, 6557, 4361, 2322, 2027, 564, and 125.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.