

Product Datasheet

Ribonuclease A (orb420130)

Description	Ribonuclease A
Conjugation	Unconjugated
Preservatives	50% (v/v) Glycerol. None
Form/Appearance	Liquid (sterile filtered)
Concentration	5 mg/mL
Storage	Store vial at 4° C or at -20° C or colder prior to opening product.
Note	For research use only
Application notes	One unit will cause an increase in the absorbance of 1.0 OD units at 260 nm at 37° C and pH 5.0 when yeast ribosomal RNA is hydrolyzed to acid soluble oligonucleotides.
Purity	Dnase - detected Protease - detected
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.