

## Product Datasheet

### Micrococcal Nuclease S7 (orb420133)

|                           |                                                                                                           |
|---------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Description</b>        | Micrococcal Nuclease S7                                                                                   |
| <b>Conjugation</b>        | Unconjugated                                                                                              |
| <b>Form/Appearance</b>    | Lyophilized                                                                                               |
| <b>Storage</b>            | Store vial at 4° C or at -20° C or colder prior to opening product.                                       |
| <b>Note</b>               | For research use only                                                                                     |
| <b>Application notes</b>  | One Unit causes a change in optical density of 1.0 at 260 nm at 37°C, pH 8.0, using DNA as the substrate. |
| <b>Purity</b>             | Purified by chromatography and assayed by SDS-PAGE.                                                       |
| <b>Hazard Information</b> | Non-Toxic                                                                                                 |
| <b>Expiration Date</b>    | 12 months from date of receipt.                                                                           |