

Product Datasheet

Histone H3 R8me2a antibody (orb420394)

Description

Histone H3 R8me2a antibody

Species/Host

Rabbit

Reactivity

C. elegans, Human, Mouse

Conjugation

Unconjugated

Tested

ChIP, DOT, IF, IHC, Multiplex Assay, WB

Applications

Immunogen

Histone H3 [Asym-dimethyl Arg8] affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic dimethylated peptide surrounding Arginine 8 of human Histone H3.2.

Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Liquid (sterile filtered)

Concentration

0.84 mg/mL

Storage

Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Histone H3 [Asym-dimethyl Arg8] antibody is tested for Western Blot, Immunocytochemistry, Immunofluorescence, Chromatin Immunoprecipitation, and Dot Blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately ~15.4 kDa corresponding to Histone H3 protein by Western Blotting in HeLa histone prep lysate or the appropriate cell lysate or extract. Epi-Plus antibody production in collaboration with Novus Biologicals.

Isotype

IgG

Clonality

Polyclonal

Purity

Anti-Histone H3 [Asym-dimethyl Arg8] was affinity purified from monospecific antiserum by immunoaffinity

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

chicken, Xenopus, Drosophila, and plant based on 100% sequence homology. Cross-