

Product Datasheet

ERK1 N-Term blocking peptide (orb420506)

Description	ERK1 N-Term blocking peptide
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Lyophilized
Concentration	1.0 mg/mL
Storage	Store vial at 2-8 °C prior to opening. Aliquot contents and freeze at -20 °C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use.
Note	For research use only
Application notes	ERK1 N-term blocking peptide is suitable as a control for polyclonal or monoclonal Anti-ERK1 in immunological assays.
Purity	Purity was determined to be greater than 95% specific peptide by analysis by HPLC and by reducing and non-reducing SDS-PAGE.
Uniprot ID	P27361
Dilution Range	WB: 50ng
Expiration Date	12 months from date of receipt.