

## Product Datasheet

### MEK1 C-Term blocking peptide (orb420509)

<b>Description</b>	MEK1 C-Term blocking peptide
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Preservatives</b>	0.01% (w/v) Sodium Azide
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	Store vial at 2-8 °C prior to opening. Aliquot contents and freeze at -20 °C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use.
<b>Note</b>	For research use only
<b>Application notes</b>	MEK1 C-term blocking peptide is suitable as a control for polyclonal or monoclonal Anti-MEK1 in immunological assays.
<b>Purity</b>	Purity was determined to be greater than 95% specific peptide by analysis by HPLC and by reducing and non-reducing SDS-PAGE.
<b>Uniprot ID</b>	<b>Q02750</b>
<b>Dilution Range</b>	WB: 50ng
<b>Expiration Date</b>	12 months from date of receipt.