

Product Datasheet

Human NONO protein (orb420632)

Description	Recombinant human Non-POU domain-containing octamer-binding protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 10xHis-taggedExpression Region: 1-471aaSequence Info: full length protein
Protein Sequence	MQSNKTFNLEKQNHTPRKHHQHHHQQQHHQQQQQPPPPPIANGQQASSQ NEGLTIDLKNFRKPGEKFTQSRSLFVGNLPPDITEEEMRKLFEKYGKAGEVFIHK DKGFGFIRLETRTLAEIAKVELDNMPLRGKQLRVRFACHSASLTVRNLPQYVSNE LLEAFSVFGQVERAVVIVDDRGRPSGKGIVEFSGKPAARKALDRCSEGSFLLTT FPRPVTVEPMDQLDDEEGLPEKLVKNQQFHKEREQPPRFAQPGSFYEYAMRW KALIEMEKQQQDQVDRNIKEAREKLEMEMEAARHEHQVMLMRQDLMRRQEEL RRMEELHNQEVQKRKQLELRQEEERRRREEEMRRQEEEMRRQEGFKGTFP DAREQEIRMGQMAMGGAMGINNRGAMPPAPVPAGTPAPPGPATMMPDGLGL TPPTTERFGQAATMEGIGAIGGTPPAFNRAAPGAEFAPNKRRRY
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	59.7 kDa
Uniprot ID	Q15233
Expiration Date	6 months from date of receipt.