

Product Datasheet

Tead1b antibody (orb420635)

Description	Rabbit polyclonal antibody to Tead1b
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	F8W5M8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42969
Uniprot ID	F8W5M8
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.