

Product Datasheet

RCC1 antibody (orb420643)

Description	Mouse monoclonal antibody to RCC1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	RCC1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	44969
Uniprot ID	P18754
NCBI	NP_001260.1 , NP_001041659.1 , NP_001041660.1 , NP_001041664.1
Dilution Range	WB: 1:2000-4000
Expiration Date	12 months from date of receipt.