

Product Datasheet

TERT antibody (orb420664)

Description	Rabbit polyclonal antibody to TERT
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	TERT
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	126997
Uniprot ID	O14746
NCBI	NP_937983.2 , NP_001180305.1
Dilution Range	IHC-P: 1:25, FACS: 1:25, WB: 1:2000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.