

Product Datasheet

RNF4 antibody (orb420671)

Description	Rabbit polyclonal antibody to RNF4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	RNF4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21319
Uniprot ID	P78317
NCBI	NP_001171939.1 , NP_002929.1 , NP_001171938.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.