



Product Datasheet

APG8B antibody (orb420675)


Description	Rabbit polyclonal antibody to APG8B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	MAP1LC3B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14688
Uniprot ID	Q9GZQ8
NCBI	NP_073729.1
Dilution Range	WB: 1:1000, IF/ICC: 1:100, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of U251 cells...



Western blot analysis of HepG2, mouse NI...



Western blot analysis of NIH/3T3, HT-108...