



Product Datasheet

CD133 antibody (orb420682)


Description	Rabbit polyclonal antibody to CD133
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PROM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	97202
Uniprot ID	O43490
NCBI	NP_001139319.1, NP_006008.1, NP_001139321.1, NP_001139322.1, NP_001139320.1, NP_001139323.1, NP_001139324.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence data is of Bone Marro...



Immunofluorescence analysis of anti-AC13...



Western blot:All lanes:Anti-hAC133-K848 ...