

Product Datasheet

HIF1Alpha antibody (orb420705)

Description	Rabbit polyclonal antibody to HIF1Alpha
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	HIF1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	92670
Uniprot ID	Q16665
NCBI	NP_001230013.1 , NP_851397.1 , NP_001521.1
Dilution Range	WB: 1:250, IHC-P: 1:50-100, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.