

Product Datasheet

GAL3ST1 antibody (orb420712)

Description	Rabbit polyclonal antibody to GAL3ST1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	GAL3ST1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48764
Uniprot ID	Q99999
NCBI	NP_004852.1
Dilution Range	WB: 1:2000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.