

Product Datasheet

SYP antibody (orb420716)

Description	Rabbit polyclonal antibody to SYP
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	SYP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33845
Uniprot ID	P08247
NCBI	NP_003170.1
Dilution Range	WB: 1:2000, IHC-P: 1:1000, IF/ICC: 1:100, FACS: 1:50
Expiration Date	12 months from date of receipt.