

Product Datasheet

DANRE pou3f3a antibody (orb420722)

Description	Rabbit polyclonal antibody to DANRE pou3f3a
Species/Host	Rabbit
Reactivity	Human, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	pou3f3a
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47216
Uniprot ID	P56224
NCBI	NP_571225.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.