

Product Datasheet

DANRE ndufaf3 antibody (orb420741)

Description	Rabbit polyclonal antibody to DANRE ndufaf3
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ndufaf3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19019
Uniprot ID	A1L1F1
NCBI	NP_001074084.1
Dilution Range	WB: 1:8000
Expiration Date	12 months from date of receipt.