

Product Datasheet

DANRE efnb2a antibody (orb420748)

Description	Rabbit polyclonal antibody to DANRE efnb2a
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	efnb2a
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36724
Uniprot ID	O73874
NCBI	NP_571098.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.