

Product Datasheet

(DANRE)edf1 antibody (orb420784)

Description	Rabbit polyclonal antibody to (DANRE)edf1
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	edf1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16075
Uniprot ID	Q6PBY3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.