

## Product Datasheet

### beta Actin antibody (HRP) (orb420797)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody against beta Actin conjugated to HRP                           |
| <b>Species/Host</b>        | Mouse                                                                                    |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | WB                                                                                       |
| <b>Target</b>              | This Beta-actin antibody is generated from a mouse immunized with a recombinant protein. |
| <b>Form/Appearance</b>     | PBS with 50% glycerol                                                                    |
| <b>Note</b>                | For research use only                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                               |
| <b>Uniprot ID</b>          | unknow                                                                                   |
| <b>Dilution Range</b>      | WB: 1:2000                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |