

Product Datasheet

EEF1A1P5 antibody (orb420809)

Description	Rabbit polyclonal antibody to EEF1A1P5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthesized Peptide
Target	EEF1A1P5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50185
Uniprot ID	Q5VTE0
Dilution Range	FACS: 1:25, WB: 1:1000-2000
Expiration Date	12 months from date of receipt.