

Product Datasheet

MYPT1 (Thr852) antibody (orb420841)

Description	Rabbit polyclonal antibody to MYPT1 (Thr852)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	PPP1R12A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	115281
Uniprot ID	O14974
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.