

Product Datasheet

AFP antibody (orb420849)

Description	Mouse monoclonal antibody to AFP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	AFP
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b,k
Clonality	Monoclonal
MW	68678
Uniprot ID	P02771
Dilution Range	IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.