

Product Datasheet

GSK3B (Phospho-Ser9) antibody (orb420856)

Description	Rabbit polyclonal antibody to GSK3B (Phospho-Ser9)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	GSK3B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46744
Uniprot ID	P49841
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.