

Product Datasheet

STREP II TAG antibody (orb420873)

Description	Mouse monoclonal antibody to STREP II TAG
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	This STREP II TAG antibody is generated from a mouse immunized with a KLH conjugated synthetic peptide of STREP II TAG.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
Uniprot ID	Unknown-STREP
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.