

Product Datasheet

NSFL1C antibody (orb420885)

Description	Rabbit polyclonal antibody to NSFL1C
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthesized Peptide
Target	NSFL1C
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40573
Uniprot ID	Q9UNZ2
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.