

Product Datasheet

BIRC6 antibody (orb420905)

Description	Rabbit polyclonal antibody to BIRC6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthesized Peptide
Target	BIRC6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	530269
Uniprot ID	Q9NR09
Dilution Range	FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.