

Product Datasheet

MIF4GD antibody (orb420909)

Description	Rabbit polyclonal antibody to MIF4GD
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthesized Peptide
Target	MIF4GD
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25423
Uniprot ID	A9UHW6
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.