

Product Datasheet

GGT3P antibody (orb420911)

Description	Rabbit polyclonal antibody to GGT3P
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthesized Peptide
Target	GGT3P
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61502
Uniprot ID	A6NGU5
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.