

Product Datasheet

LMNA antibody (orb420921)

Description	Mouse monoclonal antibody to LMNA
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	LMNA
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	74139
Uniprot ID	P02545
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.