

Product Datasheet

DLX1 antibody (orb420947)

Description	Rabbit polyclonal antibody to DLX1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthesized Peptide
Target	DLX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27320
Uniprot ID	P56177
Dilution Range	FACS: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.