

Product Datasheet

TNMD antibody (orb420955)

Description	Rabbit polyclonal antibody to TNMD
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	TNMD
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37130
Uniprot ID	Q9H2S6
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.