

Product Datasheet

alpha Synuclein antibody (orb420988)

Description	Mouse monoclonal antibody to alpha Synuclein
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	SNCA
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	14460
Uniprot ID	P37840
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.