

Product Datasheet

CCL26 antibody (orb420991)

Description	Rabbit polyclonal antibody to CCL26
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthesized Peptide
Target	CCL26
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	10648
Uniprot ID	Q9Y258
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.