

Product Datasheet

CD45 antibody (orb420999)

Description	Mouse monoclonal antibody to CD45
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant Protein
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	unknow
Dilution Range	FACS: 1:50
Expiration Date	12 months from date of receipt.