

Product Datasheet

MKS1 antibody (orb421010)

Description	Rabbit polyclonal antibody to MKS1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthesized Peptide
Target	MKS1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64528
Uniprot ID	Q9NXB0
Dilution Range	IF/ICC: 1:25, FACS: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.